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| Model Type                                | Reservoir Capacity | Temperature Range        |                         |
|-------------------------------------------|--------------------|--------------------------|-------------------------|
|                                           |                    | °C                       | °F                      |
| TC-150SD Heat only Circulator *           | 6 liters           | Ambient +10°C to 150°C** | Ambient +20° to 302°F** |
| TC-250SD Heat only Circulator *           | 10 liters          | Ambient +10°C to 150°C** | Ambient +20° to 302°F** |
| TC-550SD Refrigerating/Heating Circulator | 7 liters           | -20° to 170°C            | -4° to 338°F            |
| TC-650SD Refrigerating/Heating Circulator | 7 liters           | -20° to 170°C            | -4° to 338°F            |
| TC-950SD Refrigerating/Heating Circulator | 15 liters          | -30° to 170°C            | -22° to 338°F           |

\*, Compatible with optional TC-351 Chiller for low temperature operation to -15°C (5°F).

\*\*, Upper temperature limit based on reservoir construction. Temperature controller capable of higher temperatures.

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
|                                                            | Heating Only Circulators                  |                                     | Refrigerating / Heating Circulators |           |           |
|------------------------------------------------------------|-------------------------------------------|-------------------------------------|-------------------------------------|-----------|-----------|
|                                                            | TC-150SD                                  | TC-250SD                            | TC-550SD                            | TC-650SD  | TC-950SD  |
| Resource Disk (with Operator's Manual)                     | PS110-817                                 |                                     |                                     |           |           |
| Fitting Kit for External Applications **                   | PS 510-735 (115V)<br>PS 510-736 (230V)    |                                     |                                     |           |           |
| Bypass Tubing with Spring and Insulation                   | PS 510-711                                |                                     |                                     |           |           |
| Latex Tubing (6-ft., 0.25" ID)                             | EX-Tubing<br>(Maximum temperature 80°C)** |                                     |                                     |           |           |
| Beaker Platform(s) for Bath Reservoir<br>600 mL<br>1000 mL | PS703-038                                 | PS703-038<br>PS708-039              | PS703-038                           | PS703-038 |           |
| Deck Lid(s)                                                | PS510-781                                 | PS510-781<br>PS510-782<br>PS510-783 | PS510-784                           | PS510-784 | PS510-811 |
| RS232 Communication Output Cable                           | RSS-P1                                    |                                     |                                     |           |           |
| IEC Power Cord                                             | 1                                         |                                     |                                     |           |           |
| Certificate of Compliance                                  | 1                                         |                                     |                                     |           |           |
| Quick-Start Guide                                          | 1                                         |                                     |                                     |           |           |

\*. Includes two each: ¼" male NPT x 3/8" barb nylon fitting; ¼" male NPT x ¼" barb nylon fitting; ¼" male NPT x 3/16" barb nylon fitting; tube clamp for ¼" tubing; and tube clamp for 3/8" tubing. 230V units also include ¼" male NPT to M16 stainless steel adapter.

\*\* For higher temperatures, use part number ULA-45A

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
Serial Connector — A DB9 connector is provided on the back panel of the Controller for RS232 data communication.

| Pin | RS232 |
|-----|-------|
| 1   |       |
| 2   | TX    |
| 3   | RX    |
| 4   | DTR   |
| 5   | GND   |
| 6   | DSR   |
| 7   | RTS   |
| 8   | CTS   |
| 9   |       |

### Using Bath Controller with Rheocalc Software

The controller may communicate with Brookfield Engineering Rheocalc Software (V2.0 or higher). The controller must be connected to the computer with the appropriate RS232 serial cable.

RS232 Protocol — Rheocalc Software uses the following RS232 protocol:

Data bits — 8

Parity — None

Stop bits — 1

Flow control — None

Baud rate — Selectable (Temperature Controller/PC baud rates must match; 9600 is recommended).

### Using Bath Controller with DV-III+ or DV-III Ultra Rheometer in Standalone Mode

The bath can also be controlled directly from a DV-III+ or DV-III Ultra Rheometer (version 5.0) without using Rheocalc. Use optional cable (DVP-207) connected from the serial port on the Rheometer base to the RS232 serial port on the bath controller head. Please see additional instructions with your Rheometer Operator Manual if you wish to use this mode.

## RS232 — For use without Brookfield Software (for customers using their own software)

Serial Connector — A 9-pin D connector is provided on the back of the Temperature Controller for RS232 data communication. A serial cable that uses only the following pins should be used to connect the Temperature Controller to the computer:

- Pin #2 — Data read (data from computer)
- Pin #3 — Data transmit (data to computer)
- Pin #5 — Signal ground

**IMPORTANT:** Do not use a serial cable that connects to the unused pins on the Temperature Controller's RS232 connection.

RS232 Protocol — The Temperature Controller uses the following RS232 protocol:

- Data bits — 8
- Parity — None
- Stop bits — 1
- Flow control — None
- Baud rate — Selectable (Temperature Controller/PC baud rates must match; 57600 is recommended).

Communication Commands — Commands must be entered in the exact format shown. Do not send a LF (line feed) after the [CR] (carriage return). Be sure to follow character case exactly.

A response must be received from the Temperature Controller before another command can be sent. All responses are terminated with a single [CR], ASCII value = 13.

| COMMAND                        | FORMAT                     | RESPONSE     |
|--------------------------------|----------------------------|--------------|
| Read current fluid temperature | T{CR}                      | Tttttuz{CR}  |
| Read current set point         | S{CR}                      | Sssssuz{CR}  |
| Change set point *             | RSssssu{CR}                | RSssssuz{CR} |
| Read remaining ramp time       | M{CR}                      | Mmmmmz{CR}   |
| Set ramp time *                | RMmmmm{CR}                 | RMmmmmz{CR}  |
| Identifies controller          | I{CR}                      | Ixxxxxx{CR}  |
| Illegal command                | All illegal characters{CR} | ?{CR}        |
| Set RUN mode **                | RAz{CR}                    | RA1{CR}      |
| Set STANDBY mode **            | RAz{CR}                    | RA2{CR}      |

\*. To change set point or set ramp time, enter the command, wait for the response, and then enter RA1{CR}.

\*\* The Temperature Controller may be forced into only two states using the RA command: **RA1{CR}** (Run) and **RA2{CR}** (Standby).

Legend:

|        |                                                               |
|--------|---------------------------------------------------------------|
| ssss   | Set point (multiplied by 10)                                  |
| tttt   | Temperature (multiplied by 10)                                |
| mmmm   | Ramp time in minutes and tenths of minutes (multiplied by 10) |
| u      | Temperature units digit (F or C)                              |
| xxxxxx | Controller ID (TC-501 for all models)                         |
| z      | Controller status (1 or 2)                                    |

Controller Status [z]:

|   |                                                        |
|---|--------------------------------------------------------|
| 1 | Controller ON, in RUN mode                             |
| 2 | Controller OFF, in STANDBY mode                        |
| 3 | No probe connected to controller                       |
| 5 | Temperature reading is above or below allowable limits |
| 6 | Temperature input values outside the allowable limits  |

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## Replacement Parts. Optional Accessories & Fluids

| Replacement Parts, Standard Digital Controller                                                                | Part Number BEL                              |
|---------------------------------------------------------------------------------------------------------------|----------------------------------------------|
| Bypass Tubing Kit. Short hose with internal spring connecting the inlet and outlet ports of the circulator.   | PS510-711                                    |
| Hose Fitting Kit (SD). Includes two threaded pipe to barbed hose adapters.                                    | PS510-735 (115 volt)<br>PS510-736 (230 volt) |
| Beaker Platform, small. 1.5" height                                                                           | PS703-038                                    |
| Beaker Platform, large. 0.75"-1.5" height for use with TC-250 bath                                            | PS703-039                                    |
| Deck Lid Cover, TC-150. Flat iron opening (5.03" x 4.27")                                                     | PS510-781                                    |
| Deck Lid Covers, TC-250. Split cover male (5.69" x 5.565")                                                    | PS510-782                                    |
| Deck Lid Covers, TC-250. Split cover female (5.94" x 5.565")                                                  | PS510-783                                    |
| Deck Lid Covers, TC-250. Flat iron opening (5.03" x 4.27")                                                    | PS510-781                                    |
| Deck Lid Cover, TC-550 and TC-650                                                                             | PS510-784                                    |
| Deck Lid Cover, TC-950                                                                                        | PS510-811                                    |
| Resource Disk. Includes all literature and manuals.                                                           | PS110-817                                    |
| Reusable Air Filter, TC-550                                                                                   | PS305-057                                    |
| Reusable Air Filter, TC-650                                                                                   | PS305-054                                    |
| Reusable Air Filter, TC-950                                                                                   | PS305-055                                    |
| O-Ring for Drain Valve, TC-550 and TC-650 baths                                                               | PS400-934                                    |
| Accessory connection Hose, rated for -50° to +80°C                                                            | EX-Tubing                                    |
| Refrigeration Control Cable. TC-550, TC-650, and TC-950                                                       | PS225-651                                    |
| RS232 communication output cable                                                                              | RSS-P1                                       |
| <b>Optional Accessories</b>                                                                                   |                                              |
| Ball Valve to control external circulator flow                                                                | SB-21                                        |
| Insulated Tubing, rated for -40° to +120°C                                                                    | TC-Tubing                                    |
| High Temperature Fluoran Tubing Kit, rated for -40° to +200°C                                                 | ULA-45A                                      |
| Beaker Cover. Stainless steel split beaker cover with openings for spindle, guard leg, and temperature probe. | PB-1Y                                        |
| <b>Fluids</b>                                                                                                 |                                              |
| Algaecide, 8 oz.                                                                                              | TC-Fluid 1A                                  |
| Ethylene Glycol. Laboratory grade, 1 gallon<br>(-20° to +100°C when mixed 50/50 with distilled water)         | TC-Fluid 2                                   |
| High Temperature Fluid (50° to 150°C). 1 gallon                                                               | TC-Fluid 3                                   |
| High Temperature Fluid (100° to 200°C). 1 gallon                                                              | TC-Fluid 4                                   |
| Low Temperature Fluid (-50° to +100°C). 1 gallon                                                              | TC-Fluid 5                                   |
| Bath Cleaner, 8 oz.                                                                                           | TC-Fluid 6A                                  |